




DEFIGARD 4000

The power to save lives

With the DEFIGARD 4000, defibrillation is as easy as ABC



MULTIPULSE
BIOWAVE



SCHILLER

The Art of Diagnostics

DEFIGARD 4000

Easy and intuitive operation – one operating knob gives you access to all main functions:

- Defibrillation energy selection in manual mode, pacemaker mode selection, SAED.
- The DEFIGARD 4000 is a 200 joules biphasic defibrillator featuring *MULTIPULSE BIOWAVE*® technology.
- Selectable energy levels: 2, 4, 8, 15, 30, 50, 70, 90, 120, 150 and 200 joules.
- The bright, clear 5.9" screen is readable in all conditions and can display ECG curves, heart rate, plethysmogram and SpO₂ values.

All this is combined in a compact and light device.

Its user-friendliness, combined with a nice design and unrivalled price-performance ratio make the DEFIGARD 4000 the perfect device for all emergency situations, in the hospital environment as well as outside.



Use in the ambulance as well as on a plane:

Thanks to the mounting that is certified for 15 g, the device can be used both in the ambulance and on an aircraft.

The mounting is equipped with a DC/DC power supply, allowing the device to be charged directly via the mounting.

The device can easily be removed from the mounting.

DEFIGARD 4000 with paddles and pockets



DEFIGARD 4000

DEFIGARD 4000, the ideal defibrillator and much more:

Available with pads or paddles

- **Manual defibrillation**
- **Semi-automated defibrillation (option)**
Spoken and displayed instructions guide the user through the defibrillation process.
The defibrillation energy is automatically adjusted when children's pads are connected
- **Pacemaker (option)**
Operational mode Fix, Demand or Overdrive (frequency x3).
- **Monitoring**
1 to 6 channels - display of one ECG curve et plethysmogram OR two ECG curves (one is fix).
- **SpO₂ Masimo (option)**
Measurement is started when a sensor is detected
- **3-channel printer (option)**
- **Memory**
Up to 13 hours of ECG data or 500 events or 80 minutes of ECG data et audio recording and 500 events in SAED mode.
- **USB interface** for data transfer from the internal memory and software updates



Wall mounting





Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 709 90
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
Alfred Schiller Medical Equipment Co. Ltd.
100015 Beijing, China
Phone +86-010-52007020
diamond@schillermedical.com
www.schiller.cn



France
SCHILLER Médical S.A.S.
F-67162 Wissembourg/Cedex
Phone +33 3 88 63 36 00
Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77600 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
infoschiller@schiller-france.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85622 Feldkirchen b. München
Phone +49 89 62 99 81-0
Fax +49 89 62 99 81-54
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 (1) 383-4780 / 460-9491
Fax +36 (1) 383-4778
sales@schiller.at
www.schiller-hungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 001, India
Phone +91 22 6152 3333/ 2920 9141
Fax +91 22 2920 9142
sales@schillerindia.com
www.schillerindia.com



Japan
SCHILLER Japan, Ltd.
Hiroshima 734-8551
Phone +81 82 250 2055
Fax +81 82 253 1713
koji.maekawa@schiller.jp
www.schiller.jp



Croatia
Schiller medicinski instrumenti d.o.o.
10000 Zagreb
Phone +385 1 309 66 59
Fax +385 1 309 66 60
info@schillerzg.hr
www.schiller.ch



Latin America
SCHILLER Latin America, Inc.
Doral, Florida 33172
Phone +1 954 673 0358
Fax +1 786 845 06 02
info@schillerla.ch
www.schillerla.ch



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 8432089
Fax +48 22 8432089
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
SCHILLER AG Rep. office
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 970 11 33
mail@schiller-ag.com
www.schiller-cis.com



Serbia
SCHILLER d.o.o.
11010 Beograd
Phone +381 11 39 79 508
Fax +381 11 39 79 518
info@schiller.rs
www.schiller.rs



Slovenia
SCHILLER d.o.o.
2310 Slovenska Bistrica
Phone +386 2 843 00 56
Fax +386 2 843 00 57
info@schiller.si
www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
E-28230-Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es




Switzerland
SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
Okmeydanı-Sisli - Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com



USA
SCHILLER America Inc.
Doral, Florida 33172
Phone +1 786 845 0620
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com



SCHILLER

The Art of Diagnostics